



BACKGROUND

- ⁶⁸Ga-DOTATATE PET/CT is a promising imaging tool used to detect and monitor disease in patients with advanced gastroenteropancreatic neuroendocrine tumors (GEP-NETs). (Figure 2, examples of intersection method Figure 3) Patients often present with a high disease burden, sometimes with tens to hundreds of lesions, making comprehensive lesion-wise assessment clinically infeasible. 8.0 <u>≦</u> 25 -Here, we implement automated, convolutional neural networkbased (CNN) methods for automatic individual lesion \$ 20 · detection and disease burden assessment. MATERIALS AND METHODS ษี 0.4 Baseline & follow-up ⁶⁸Ga-DOTATATE PET/CT images from 59 patients with **⊂** 0.3 GEP-NETs undergoing theranostic ¹⁷⁷Lu-DOTATATE (Lutathera) therapy were S 0.2 retrospectively analyzed (116 total scans, 1-7 per patient). 0.1 • Two different CNNs (Figure 1) were trained separately using 5-fold cross validation • Ensembles of the two CNNS were created with the intersection and union of the outputs Figure 2: Detection performance of all methods **(A)** • Patient-level SUV_{max} was correctly captured in 49 of 59 scans 128x28x28x48 256x14x14x24 Patient-level SUVmax Patient-level SUVmean 256x14x14x24 Semantic Seq. Features cl Classification Feature Pyramid Network = 3x3x3 conv – Instance Normalization – Leaky ReLu = Input = 1x1x1 conv - softmax Image from Jaeger et al., 2018 5 150 · Figure 1: Architecture for the (A) Retina U-net [1] and (B) nn-Unet [2] model used in this study. **Performance Quantification** 50 - Lesion detection performance quantified with lesion detection sensitivity and 49/59 pts with correct SUVmax number of false positives (FPs) per patient Quantification performance of SUV_{max}, SUV_{mean}, SUV_{total}, and total volume Expert-based Expert-based assessed with Pearson's correlation coefficient (R).





Automated detection and quantification of neuroendocrine tumors on ⁶⁸Ga-DOTATATE PET/CT images using a U-net ensemble method

Amy J Weisman¹, Ojaswita Lokre¹, Brayden Schott², Victor Santoro-Fernandes², Robert Jeraj^{1, 2}, Timothy G Perk¹, Steve Y Cho², Scott B Perlman² ¹AIQ Solutions, ²University of Wisconsin-Madison

RESULTS

• A total of 2,634 lesions from the 59 baseline PET/CT images were contoured by a radiographer (range: 1-239 lesions per scan). • Median sensitivity was 87% with 2 FPs/patient for nnU-net, and 92% sensitivity with 5 FPs/patient for retina U-net (Figure 2) • The union ensemble achieved 93% sensitivity with 5 FPs/patient, and intersection achieved 82% sensitivity with 2 FPs/patient






Figure 3: Intersection method performance in 3 patients

----- True Positives ------ False Positives ------ False Negatives

• For baseline patient-level quantification, the ensemble intersection method achieved the best overall quantification performance, with Pearson correlation coefficients of R=0.95 for SUV_{mean}, R=0.97 for SUV_{total}, and R=0.92 for total volume (**Figure 4**).




Figure 4: Patient-level quantification performance for the intersection method



KEY FINDINGS

An ensemble of two U-net based CNNs trained for lesion detection achieved excellent performance for quantifying patient-level PET imaging metrics.

Despite a lower sensitivity, the method with the fewest false positives achieved the best quantification performance, indicating the majority of missed lesions have low uptake and represent a small fraction of the total disease burden.

REFERENCES

[1] Jaeger, P. F. et al. Retina U-Net: Embarrassingly Simple Exploitation of Segmentation Supervision for Medical Object Detection. Proc. Mach. Learn. Res. NeurIPS 2019 1–12 (2018).

[2] Isensee, F. *et al.* nnU-Net: a self-configuring method for deep learning-based biomedical image segmentation. Nat. *Methods* **18**, (2020).

DISCLOSURES

Authors AJW, OJ, and TGP are employed by AIQ Solutions. RJ is a cofounder and board member of AIQ Solutions. SBP is a consultant for AIQ Solutions. Remaining authors have no relevant relationships to disclose.

CONTACT INFORMATION

Amy.Weisman@aiq-solutions.com



